Lumity Housing Colours



Instant Fit

Domes







Vented Dome

*Instant fit domes are not compatible with MAV Receivers. Universal SlimTips are available as an instant fit or trial solution for MAV receivers.

Custom Ear Pieces

























(solid hard)








SlimTip Titanium AV (MAV receivers only)

ActiveVent Receiver (MAV)

The world's first intelligent hearing aid receiver


- Intelligently balances speech clarity in noise and listening comfort
- Provides superior sound quality for streamed media
- Compatible with all Lumity rechargeable hearing aids
- Compatible with universal tips and SlimTip Titanium







Acoustic Coupling




M Receiver Gain 51(2cc) • MPO 114(2cc)


MAV Receiver Gain 50(2cc) • MPO 119(2cc)



125 250 500 1k 2k 4k 8k Mild to moderately-severe hearing loss

P Receiver Gain 59(2cc) • MPO 122(2cc)






Moderately severe to profound hearing loss

More Phonak technology to enhance your experience







Marvel



Paradise







CROS Paradise



Behind-the-ear



312




312 NW 0



10 NW 0







Phonak Serenity Choice Hearing Protection

- Universal and custom-made hearing protection designed to cancel noise and loud sound
- Keep relevant sound and speech
- Lets the ear breathe



sonova
















PHONAK life is on

Features at a Glance

Phonak Lumity.

Phonak Lumity Models

Audéo™ L

A Sonova brand



Audéo L-R

Rechargeable



Rechargeable Charging Options



Audéo L-RL P90/P70/P50/P30 Rechargeable Waterproof*



Adjust volume and programs with the easy to use remote control



Wireless Accessories


Phonak PartnerMic™ Lapel style directional microphone for better one to one conversations





headphones







Charger Case Ease

Contact charging

Magnetic retention

and Audeo L-RT

• 3 hours to fully charge

Compatible with Audeo L-R

Roger X Receiver Wirelessly install Roger X into Lumity hearing aids for RogerDirect




Table Top, Lanyard,

& Pointing mode

Charger Case Go

Induction charging

Compatible with Audeo L-RL

Built-in battery for 3 full

charges on the go



Roger Select Table Top, Lanyard,

& Select mode

Roger TableMic II Table Top mode




The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. 028-2146-02/2022-08/ Printed in Canada © Phonak Canada All rights reserved

Introducing SmartSpeech Technology

A collection of features working together to focus on improving speech understanding in all envrionments:

- Understand more speech in quiet and noise
- Spend less effort on listening





The first fully rechargeable product range from Phonak - the future is on

Better for you

- Access additional hearing performance features only available because of rechargeable technology
- Easy to use

Better for the environment

• Save on average 800 disposable batteries during the lifetime of a pair of hearing aids

vels



BENEFITS

Provides significant improvement in speech understanding in noise and far

Seamlessly adjusts for more speech focus in the front in challenging noise

Intelligently balances hearing performance and comfort of own voice

In noise, zooms in on the voices you are facing

Automatically adjusts to optimize hearing while walking

Boosts soft and distant speech in quiet environments

Provides listening comfort and clarity in quiet

Adapts to optimize speech in noisy environments

Reduces car and road noise

Improves hearing from the side and behind when in noise

Listening comfort in loud situations when no speech is present

Uses StereoZoom 2.0 to help focus on speech in a loud noise environment

Offers a full and rich sound quality to your favorite tunes

Reduces disturbing echo so you can hear speech clearly

Suppresses wind noise so you can enjoy time outdoors

Reduces background noise without reducing speech

Cushions sudden loud or unpleasant impulse sounds

Removes echo/reverberation to provide a cleaner sound

Noise generator to reduce the impact of your tinnitus


Pause/resume streaming, accept/end calls, and access voice assistant

Tracks and focuses on the dominant speaker from different sides

Eliminates whistling from the hearing aids before it is audible

Reduces background noise around you for improved hearing









FEATUREC

For feature details and power levels per model, please consult Target Software







Fully participate in challenging group conversation	ons
Enjoy hearing aids that automatically adjust	

	Dynamic Noise Cancellation	NEW in 70 tech		Speech Enhancer
ActiveVent Receivers		0 0	0	((o)) UltraZoom
Motion Sensor Hearing	((_k))	00	0	NEW StereoZoom 2.0
	Roger Direct		() (()	SpeechSensor

Better for your loved ones

• Safer use - no batteries for small children or pets to accidently get hold of

Premium

eatures that shine	in	all	technol	ogy	leve
				٠,	

FEATURES

ActiveVent Performance

Motion Sensor Hearing

Dynamic Noise Cancellation

Roger Technology

Speech Enhancer

SpeechSensor

Calm situation

Speech in Noise

Comfort in Noise

Speech in Loud Noise

Speech in Car

WindBlock

WhistleBlock

NoiseBlock

SoundRelax

EchoBlock

Tinnitus Balance

Tap Control

Comfort in Echo

Speech in 360° (manual)

Music

UltraZoom

™ ®	Universal
$\mathbf{X} \square$	Bluetooth
	connectivity

DENIELITO

	martSpeech
dard	S 5.0 S
ano	AutoSense 09
7	

	FEATURES	BENEFITS
	ActiveVent Performance	Intelligently balances hearing performance and comfort of own voice
SmartSpeech	Roger Technology	Provides significant improvement in speech understanding in noise and far distances
mart	UltraZoom	In noise, zooms in on the voices you are facing
S	Motion Sensor Hearing	Automatically adjusts to optimize hearing while walking
S 5.0	Calm situation	Provides listening comfort and clarity in quiet
AutoSense OS	Speech in Noise	Adapts to optimize speech in noisy environments
AutoS	Comfort in Noise	Listening comfort in loud situations when no speech is present
	Music	Offers a full and rich sound quality to your favorite tunes
	WindBlock NEW in 50 tech	Suppresses wind noise so you can enjoy time outdoors
	WhistleBlock	Eliminates whistling from the hearing aids before it is audible
	NoiseBlock	Reduces background noise without reducing speech
	Tinnitus Balance	Noise generator to reduce the impact of your tinnitus

0 d Van 0

	FEATURES	BENEFITS
	ActiveVent Performance	Intelligently balances hearing performance and comfort of own voice
Autosense Os 5.0 Smartspeech	Roger Technology	Provides significant improvement in speech understanding in noise and far distances
	UltraZoom	In noise, zooms in on the voices you are facing
	Motion Sensor Hearing	Automatically adjusts to optimize hearing while walking
	Dynamic Noise NEW in 70 tech	Reduces background noise around you for improved hearing
	StereoZoom 2.0 NEW	Seamlessly adjusts for more speech focus in the front in challenging nois
	Calm situation	Provides listening comfort and clarity in quiet
	Speech in Noise	Adapts to optimize speech in noisy environments
	Comfort in Noise	Listening comfort in loud situations when no speech is present
	Music	Offers a full and rich sound quality to your favorite tunes
	Speech in 360° (manual)	Tracks and focuses on the dominant speaker from different sides
	Speech in Loud Noise (manual)	Uses StereoZoom 2.0 to help focus on speech in a loud noise environmen
	WindBlock	Suppresses wind noise so you can enjoy time outdoors
	WhistleBlock	Eliminates whistling from the hearing aids before it is audible
	NoiseBlock	Reduces background noise without reducing speech
	SoundRelax	Cushions sudden loud or unpleasant impulse sounds
	Tinnitus Balance	Noise generator to reduce the impact of your tinnitus
	Tap Control	Pause/resume streaming, accept/end calls, and access voice assistant

For feature details and power levels per model, please consult Target Software For feature details and power levels per model, please consult Target Software