

ReSound OMNIA · Made for hearing in noise

ReSound GN

Want your patients to hear
their best in noise and enjoy
sounds all around them?
Now they can.



Your patients can now hear their
best in noise with an incredible 150%
improvement in speech understanding.*
Powered by the industry's only
Organic Hearing™.




GN Making Life Sound Better

*4.3dB SNR improvement to speech understanding compared to legacy products. Jespersen et al (2022)




Wow your patients from the start

Hearing in noise is a huge challenge for people with hearing loss. We have taken a great leap forward with the sound benefits that ReSound OMNIA offers, which help people hear what they want even better in noise. You now have everything you need to impress your patients with ease.



Powered by Organic Hearing

ReSound OMNIA is powered by the industry's only Organic Hearing, our unique approach to creating hearing solutions that sound natural, feel natural and connect your patients naturally to the world.

-  Sounds natural
-  Feels natural
-  Connects naturally



How did we achieve it?

Combining the strength of humans and machine with our cutting-edge technology and leading audiological expertise, we developed our most intelligent and automated steering and beamforming system yet. This gives users superior hearing clarity and understanding in all sound environments, so they can love the everyday details in life again.

1 

Widening our sound environment mapping to 360 All Around

2 

Speeding up our processing and steering for faster mode detection and selection

3 

Narrowing the beamforming to give even more precision for our best speech clarity in any environment


It all happens simultaneously

All working at the same time across any sound environment, enabling your patients to automatically tune into what they want to listen to and not have the hearing aids choose for them.

The everyday details that make a huge difference

ReSound OMNIA is made for hearing in noise and fits naturally into people's lives, starting with the popular Receiver-in-Ear (RIE) style.

- Our best speech understanding in noise*¹
- Hearing clarity all around¹
- Natural own voice perception²
- Optimal comfort on the ear
- All-day power with up to 30 hours use on one charge
- All-weatherproof design for less maintenance
- Check My Fit to help patients ensure their hearing aids sit correctly on the ear
- Crystal clear phone calls
- Hands-free calls for iPhone and iPad**
- Direct streaming from your favorite devices



* Compared to our legacy technology
¹ Jespersen et al (2022)
² Groth et al (2022)

Learn more at pro.resound.com

Improve onboarding and usage

The new and unique Check My Fit feature helps users to wear their hearing aids properly and get the most out of ReSound OMNIA each day, which can lead to fewer service visits.

Available in the ReSound Smart 3D™ app.




** ReSound OMNIA hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), and iPad mini (6th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.




Key feature overview

	Feature Name	Tech level 9	Tech level 7	Tech level 5
Comfort	Noise Tracker II	5 settings	3 settings	2 settings
	Wind Guard	3 settings	2 settings	On/Off
	Impulse Noise Reduction	3 settings	3 settings	On/Off
	Expansion	3 settings	2 settings	On/Off
	Environmental Optimizer II	•		
	Environmental Optimizer		•	•
Speech Understanding	360 All-Around	•		
	All Access Directionality		•	
	Binaural Directionality III			•
	Front Focus	•		
	Ultra Focus		•	
	AutoScope Adaptive Directionality	•		
	MultiScope Adaptive Directionality		•	
	Synchronized Soft Switching	•	•	•
	Adaptive Directionality	•	•	
	Fixed Directionality	•	•	•
	M&RIE	•	•	•
	Spatial Sense	•	•	•
	Adjustable Directional Mix	•		
Feedback Management	DFS Ultra III	•	•	•
	Music Mode	•	•	•
	Auto DFS	•	•	•
Convenience	Ear-to-ear synchronization	•	•	•
	Hands-free calls for iPhone & iPad	•	•	•
	Made for iPhone, iPad and iPod touch (MFi)	•	•	•
	Android streaming (ASHA)	•	•	•
	Wireless Accessories	•	•	•
Fitting Flexibility	Max Number of Programs	4	4	4
	ReSound Assist	•	•	•
	ReSound Assist Live	•	•	•
	Low Frequency Boost	3 settings	2 settings	On/Off
Amplification strategy (WDRC-Semi-Linear-Linear)	3 settings	3 settings	2 settings	



Stylish colors

Available in eight beautiful tones for you to help guide your patients in choosing a color to match their preferences.



Sparkling Silver Champagne



Gold Bronze



Warm Grey Graphite



Espresso Deep Black

Introducing improved and more efficient ReSound Smart Fit™

Provide your best hearing care, even remotely, fitting many needs and preferences efficiently with our optimized navigation flow and fitting process in ReSound Smart Fit 1.15.

Manufacturer according to FDA:

GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-248-4327
pro.resound.com

ReSound Government Services

8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
gs.resound.com

Manufacturer according to Health Canada:

ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com

© 2022 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple and the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC

